

CLIFTON NE1V OVEN/ INCUBATOR

Model: 800-NE1V.32E-00/ 800-NE1V.56E-00



■ **If you're looking for a dependable oven/incubator with a controlled temperature environment for various applications, the NE1V is an excellent choice.**

With its digital PID controller, temperature regulation is precise and reliable. Plus, the automatic alarm feature adds an extra level of safety. And, with its insulated chamber, energy consumption is kept low.

- Dual temperature display actual and set values
- Stainless steel chamber
- Gravity convection chamber
- Exhaust outlet
- Lockable and insulated door
- Silicone door chamber seal
- Anti-bacterial paint exterior
- Levelling feet, adjustable in each corner
- Factory two-point calibration, 37°C and 150°C
- Two-year warranty
- Made in Britain



G A L L E N K A M P

SPECIFICATIONS

NEIV-56

CAT. No.	800-NEIV.32E-00	800-NEIV.56E-00
Temp. Range	Ambient +5°C -250°C	
Temp. Control	PID Digital, LED Display Showing Actual & Set Temperature, Resolution 0.1°C	
Temp. Sensor	PT100 With Diagnostics Check On Power Up	
Temp. Fluctuation 37C	±0.3°C	
Temp. Variation 37C	±1.0°C	
Temp. Fluctuation 150C	±0.5°C	
Temp. Variation 150C	±2.0°C	
Temp. Heat Up 150C	30 Minutes	40 Minutes
Temp. Control With Timer Program	Run Back Timer 00:01 to 99:59 (HH:MM) At Set Temperature	
Temp. Alarm	Over Temperature Alarm Automatically Set 5°C Above Set Temperature	
Temp. Door Open Recovery 150C	16 Minutes	
Over Temp. Limiter	Class 1	
Shelf no. Supplier/ Max	1 Wire / 2 Max	2 Wire / 4 Max
Shelving Locations	Available 80mm Spacing	
Heater Power	1050W	
Voltage	230V 50/60Hz	
Energy Use At 37C	0.02kw/h	
Energy Use At 150C	0.23kw/h	

ADDITIONAL INFORMATION

CAT. No.	800-NEIV.32E-00	800-NEIV.56E-00
(Internal) WxDxH MM	400x270x320	400x390x360
(External) WxDxH MM	740x460x560	740x610x600
Chamber Capacity	32L	56L

ACCESSORIES

Description	CAT. No.
Wire shelf stainless steel for 32L	WS-32
Wire shelf stainless steel for 56L	WS-56

ORDERING INFORMATION

Description	CAT. No.
Clifton NEIV Oven/Incubator 32 lts 230V	800-NEIV.32E-00
Clifton NEIV Oven/Incubator 56 lts 230V	800-NEIV.56E-00



G A L L E N K A M P

Beacon House
Nuffield Road
Cambridge
CB4 1TF
United Kingdom

Tel: +44(0)1223 727 100
Email: sales@fistreem.co.uk
Web: www.gallenkamp-labs.com